




Designed for Revisions. Designed for Surgeons.

Visit our website at [www.exodusrevision.com](http://www.exodusrevision.com) for more information

00-507-400-18/19	1	Exodus Handle - Strikeplate
00-507-400-16	1	Exodus Handle
00-507-400-40	1	Tommy Bar
00-507-400-76	1	Jack Rabbit Slap Hammer
00-507-400-77	1	Jack Rabbit Slap Hammer Sling
00-507-450-00	2	Ball Nose T40 Hex Driver
00-507-401-00	2	Cup Cutter
00-507-401-28	1	Cup Cutter Ball 28
00-507-401-32	1	Cup Cutter Ball 32
00-507-401-36	1	Cup Cutter Ball 36
00-507-401-40	1	Cup Cutter Ball 40
00-507-401-44	1	Cup Cutter Ball 44
00-507-430-41	1	Drill Centering Spacer 10mm
00-507-430-42	1	Drill Centering Spacer 12mm
00-507-430-43	1	Drill Centering Spacer 16mm
00-507-401-90	2	Hudson Bell Blocker, Orange

### Handle Tray:

- The Exodus Revision Hip System is designed for **minimal** bone loss with **minimal** effort and **minimal** damage to the surrounding tissues.
- Modular handles to interface with full suite of attachments.
- The Exodus Strikeplate handle has a curved strike-plate undersurface that enables a firm grip between the palm and fingertips while back-striking during extraction. The silicone grip on the handle is designed to reduce fatigue while providing grip, ergonomics and comfort during use.
- When significant impact loads are required, the Exodus slap hammer is uniquely designed to controlled extraction forces.




## HANDLE TRAY

00-507-481-01	*	Medial Stem Osteotome, SU
00-507-481-02	*	Lateral Stem Osteotome, SU
00-507-481-03	*	Chisel, Long, SU
00-510-482-01	2	SM Cup [48-58]-Short Talon, SU
00-510-482-02	2	SM Cup [48-58]-Long Talon, SU
00-510-483-01	2	LG Cup [60-74]-Short Talon, SU
00-510-483-02	2	LG Cup [60-74]-Long Talon, SU
00-510-488-01	2	Threaded Stem Extractor #1, SU
00-510-488-02	2	Threaded Stem Extractor #2, SU
00-510-488-03	2	Threaded Stem Extractor #3, SU
00-510-488-04	2	Threaded Stem Extractor #4, SU
00-510-484-01	1	Goat Foot, SU
00-510-484-03	1	Chisel, Short, SU
00-510-484-07	1	Chisel, Medium, SU
00-510-485-04	2	Trunnion Locker, SU
00-510-486-01	2	Knee Claw, SU
00-510-487-01	1	Straight Gouge, SU
00-510-487-02	1	V-Splitter, SU
00-510-487-03	1	Reverse Curette, SU
00-510-487-04	1	Tap, SU
00-510-487-05	1	Cement Drill, SU
00-510-489-01	1	Short Hockey Stick Right, SU
00-510-489-03	1	Long Hockey Stick Right, SU

### Single Use (SU) Attachments:

- Single-use instruments ensure that the attachments are sharp and free of damage
- The instruments provide options to both break cement interfaces and cut through in-growth along the implant contour for easy removal.

*\* Note these items are sold in a Single Use Femoral Hip Stem Removal Sterile Pack: 00-510-400-00 and distributed through [ZIMMER BIOMET](http://www.zimmerbiomet.com) directly.*



## ATTACHMENT TRAY